



APPLICATION TECHNIQUE STERILIZATION


Mr. Edmond Ng

Sales & Marketing Director – IdsMED Việt Nam

Demands on the Sterilization process “ PAST”




- Technical construction according EN 285 (1996) „Type proof”
- Demands of the customer for the sterilization process
- Reached limits of the sterilizer



Demands on the Sterilization process “NOW”



- Technical construction according EN 285 (2006) „Type proof”
- Reached limits of the sterilizer
- Demands of the customer for the sterilization process „all-in-one device suitable for every purpose” and this for no costs“




Customer requirements



- Sterile medical products
- Short batch times
- Maximum load
- No residual moisture
- All container and packaging possibilities are unlimited usable
- Everything has to be possible

Factors for a good Sterilization process



Packaging



Rigid Packaging

- Container built according to EN 868 Part 8
- Test load 10 kg per StU

Drying performance according to EN 868 Part 6

- Residual moisture < 0.2% (Metal load)
- Residual moisture < 1.0% (Textile load)



Material properties



- Container and cover made of stainless steel
- Container and cover made of Aluminium / Metal trays
- Aluminium container with stainless steel cover / Metal trays
- Aluminium container with stainless steel cover / Metal trays / Condensate drain
- Aluminium container with plastic cover / Metal trays
- Container and cover made of Aluminium / Plastic trays
- Aluminium container with plastic cover / Plastic trays

Calculation Amount of condensate (Stainless steel)



Heat requirement ($Q = m \times c \times \Delta t$)

$$= 12\text{kg} \times 0.5\text{kJ/kg} / \text{K} \times 114\text{K} = 686 \text{ kJ}$$

Δt Temperature difference ($20^\circ\text{C} - 134^\circ\text{C}$) c Specific heat capacity (1kJ for 1°C)

Energy capacity of steam 1 kg Steam \cong 2000kJ (r)

Steam quantity (md) = Quantity of condensate $Q = md \times r$

$$md = Q/r = 686\text{kJ} / 2000\text{kJ} = 0.342 \text{ kg} = 343\text{g Steam}$$

Calculation Amount of condensate (Plastic)



Heat requirement ($Q = m \times c \times \Delta t$)

$$= 12\text{kg} \times 2.5\text{kJ/kg} / \text{K} \times 114\text{K} = 3420\text{kJ}$$

Δt Temperature difference ($20^\circ\text{C} - 134^\circ\text{C}$)

Energy capacity of steam 1 kg Steam \cong 2000kJ (r)

Steam quantity (md) = Quantity of condensate $Q = md \times r$

$$md = Q/r = 3420\text{kJ} / 2000\text{kJ} = 1.71\text{kg} = 1710\text{g Steam}$$

**Plastic needs
approx. 5 x more energy to
heat up!**

- More Energy**
- = More condensate**
- = Longer drying time**
- = Longer process time**




Flexible Packaging

- Papers (plain, crepe, semi crepe, fleece)
- Foils / pouches (Gas, Steam, Gamma-rays)
- No textiles (Textile alternative)
- Textiles



Materials Characters

- Cotton
- Cellulose 100%
- Cellulose with binding material
- 100% Polypropylene
- Polypropylene with absorber coating




Norm: EN 285 (Residual moisture)

- Sterilization process 3.5min at 134°C
- Drying capacity according to EN 285 (8.4.2)
- Test container 1 StU (14.1 kg ± 0.4 kg)
(Container (4,2 kg) + Screws (8,6 kg)+ Tray (1,3 kg) + Cotton wrapping)
- Residual moisture < 0.2% (28.2 g with metal load)



Condensate within the wrapping

Increase of weight in %




If the weight of the container increases, also the amount of the condensate increases!

Reminder

With increase of the container weight the condensate amount increases!

More container weight
= More condensate
= Longer drying time
= Longer process time



Condensate residue (after drying)

Visible

- Drops
- Water accumulation




Non visible

- Integrated in sterilisation fleece




Example with 10 ml water

- 100% Cellulose with chemical binding material



- Textile



Example with 10 ml water

- Polypropylene with absorber




- 100% Polypropylene Cellulose



Visual check (Dryness)

- Drop < 1ml
- Big drop < 5ml
- Accumulation approx. 10ml
- Allowed remaining quantity according to Norm: **28g**

28g = WET !!



Picture source of error



Picture Container full



- No maximum load indication from container producer



Picture Container condensate

- After a load of 24kg instruments
- Marking from filter lid



How to do it right: Correct loading



- Light sterile goods in flexible packs should be loaded on the top.
- Heavy containers should be placed at the bottom of the chamber.

Reasons:

- The heavier the items are, the higher is the amount of condensing water
- Lighter packs must be “protected” against drops of condensing water



How to do it right: Correct loading

- Pouches



- Unhindered in – and out of air and steam

Bad examples

- Overloaded, packaging can not breath



- Foils down side



Bad examples

- Textiles in contact with sterilizer chamber
(Loading index = width of batch cart)



Influence of packing materials to drying result

Packing material	Drying result	Type
Cotton	++	1- layer
Cellulose	++	1- layer
Cellulose	+	2- layer
Polypropylene + mixed fibres	++	1- layer
Polypropylene + mixed fibres	+	2- layer
100% Polypropylene	O	1- layer
100% Polypropylene	-	2- layer

--	-	O	+	++
very bad	bad	average	good	Very good

Different containers led to different drying results



Container	Drying result
Container and cover made of Aluminium / Metal trays	++
Aluminium container and cover made of stainless steel / Metal trays	+
Aluminium container with plastic cover / Metal trays	O
Container and cover made of Aluminium / Plastic trays	-
Aluminium container with plastic cover / Plastic trays	--

--	-	O	+	++
very bad	bad	average	good	Very good

Residual moisture MATRIX Container-10kg



▪Based on Sterilization process 3.5min. 134°C, 25 Min drying

Packaging / Container	Cotton cloth	100 % Cellulose	Polypropylene and mixed fibres	100% Polypropylene
Sainless steel with cond. drain	😊	😊	😊	😊
Sainless steel w/o cond. drain	😊	😊	😊	😊
Aluminium with Aluminium cover	😊	😊	😊	😊
Aluminium with plastic cover	😊	😊	😊	😊

bad	average	good	Very good
		😊	😊


Residual moisture MATRIX Container 12kg - 14kg



Packaging	Cotton cloth	100 % Cellulose	Polypropylene and mixed fibres	100% Polypropylene
Container				
Sainless steel with cond. drain	😊	😊	😊	
Sainless steel w/o cond. drain	😊	😊		
Aluminium with Aluminium cover	😊	😊		
Aluminium with plastic cover	😊	😊		
	bad	average	good	Very good
			😊	😊

Process


No.	CYCLE	STERILIZE TEMP	STERILIZE TIME	DRY TIME	RECOMMENDED LOAD
1	Instruments 134°C	134°C	4 minutes	25 minutes	Instrument intrays or containers max. 7-8kg/ tray or container



Process



No.	CYCLE	STERILIZE TEMP	STERILIZE TIME	DRY TIME	RECOMMENDED LOAD
7	Heavy Instruments 134°C	134°C	4 minutes	20 min. + 3*1,5 min	Heavy instruments in stainless steel containers, max. 14kg / container (gross weight)




Optimizing Residual moisture

- Total weight of container shall not exceed. 14 kg per StU
- Filling level of container max. 3cm below upper edge
- Minimize the use of plastic trays
- Use well absorbing wrapping material
- Pulsed-vacuum-dry is more effective vacuum dry



Source of errors: Loading

- Container is too small
- Filling level of container is exceeded
- Size of packaging paper is too big
- Stacked containers
- Horizontal positioning of flexible pouches



Source of errors: Installation

Steam to wet

- Steam supply without condensate trap
- Power of steam generator not sufficient
- Condensate trap (mechanical better than thermal)





Xin cảm ƠN

Mr. Edmond



We care with a
HEART
because your life matters