
3M Science.
Applied to Life.™

3MSM Health Care Academy

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
NHIỄM KHUẨN HUYẾT QUA CATHETER (CRBSI = Catheter-related bloodstream infection)

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

- Mermel LA. Prevention of intravascular catheter-related infections. *Ann Intern Med* 2000;132:391—402.
- CDC. National Nosocomial Infections Surveillance (NNIS) System report, data summary from October 1986--April 1998, issued June 1998. *Am J Infect Control* 1998;26:522--33
- Digiovine; The attributable mortality and costs of primary nosocomial bloodstream infections in the intensive care unit. *Am J Respir Crit Care Med* 1999;160:976—81.
- Digiovine B, Chenoweth C, Watts C, Higgins M. The attributable mortality and costs of primary nosocomial bloodstream infections in the intensive care unit. *Am J Respir Crit Care Med* 1999;160:976—81

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Organisms at the insertion site (extraluminal) are source of infection in the majority of cases in 'early' infections

Mermel L, Clin Infect Dis (2011) 52: 211-212

Recommendations, Guidelines & Standards

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Centers for Disease Control and Prevention (CDC)

- 2011 Guidelines for the Prevention of Catheter-Related Bloodstream Infections

Society of Infusion Nurses (INS)


- 2016 Standards of Practice for Infusion Nurses

Society for Healthcare Epidemiology of America (SHEA)

- 2014 Strategies for the Prevention of Bloodstream Infections Associated with Central Lines in Acute Care Hospitals

The Association for Professionals in Infection Control and Epidemiology (APIC)

- 2016 Guide for prevention of Central Line Associated Bloodstream Infections (CLABSI)

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
Chlorhexidine-impregnated dressings with an FDA-cleared label that specifies a clinical indication for reducing catheter-related bloodstream infection (CRBSI) or catheter-associated blood stream infection (CABSI) are recommended to protect the insertion site of short-term, non tunneled central venous catheters. (Category IA)8-12

1. Use either sterile gauze or sterile, transparent, semipermeable dressing to cover the catheter site. (IA) – page 42
2. Use a chlorhexidine-impregnated sponge dressing for temporary short-term catheters in patients older than 2 months of age if the CLABSI rate is not decreasing despite adherence to basic prevention measures, including education and training. (IB) - page 44

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
Align your protocols with standards of practice.

The Centers for Disease Control and Prevention (CDC), Infusion Nurses Society (INS), Association for Professionals in Infection Control and Epidemiology (APIC), Society for Healthcare Epidemiology of America (SHEA), and other organizations offer evidence-based best practices to help minimize I.V. site complications. Choose 3M™ Tegaderm™ CHG Dressings and be sure you're meeting or exceeding best practices for better patient and economic outcomes.

Site visibility

The CDC and *Infusion Therapy Standards of Practice* recommend the use of transparent dressings because they permit continuous visual inspection of the catheter site.^{3,4}

Antimicrobial protection

APIC, INS and the SHEA Compendium recommend CHG-impregnated dressings.^{4,5,6} In use for over 50 years, CHG has proven to be an effective antimicrobial. Bacterial resistance to CHG has been rare.⁷

Consistent application

The International Organization of Standards promote the importance of medical device design to support correct use, patient safety, user satisfaction and to reduce medical device-related errors.⁸


Catheter securement

The CDC and *Infusion Therapy Standards of Practice* recommend the use of sutureless securement devices to minimize the risks of movement, dislodgement, and needlestick injuries.^{3,4}

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FIGURE 7: MOST COMMON PATHOGENS CAUSING CLA-BSIs⁵

Tegaderm™ CHG Dressing demonstrates *in vitro* efficacy against

Enterococcus
(5 strains)

Pseudomonas
aeruginosa
(5 strains)

Candida
(2 strains)

Staphylococcus
aureus (8 strains)

Escherichia
coli (1 strain)

Coag Neg
Staph (7 strains)

Klebsiella
(2 strains)

Enterobacter
(1 strain)


Other
(6 strains)

37 strains of microorganisms

including gram positive bacteria, gram negative bacteria and yeast.

Figure 1: Catheter-related bloodstream infections


Hensler J (2009), et al. Growth inhibition of microorganisms involved in CA-infections by an antimicrobial transparent IV dressing containing CHG. European Society of Clinical Microbiology and Infectious Diseases, May 2009.

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Application

1. Allow all antiseptics and skin protectants to dry completely.

2. Peel the liner from the dressing. Turn over so the adhesive faces the skin. Ensure that the CHG gel pad covers the catheter insertion site (and suture sites when possible).

3. Be careful not to stretch the dressing at placement. Apply firm pressure to the entire dressing starting over the gel pad to enhance adhesion.

4. Apply pressure to securement border with one hand, while removing paper frame with opposite hand.

5. Remove notched tape strip from paper frame. Orient the notch towards the dressing and apply the tape strip under the extension tubing and over the dressing border. Remove adhesive-free tabs.

6. Document the dressing change information on the label strip. Apply label strip on top of dressing, over catheter lumen(s). Remove adhesive-free tabs.

Monitoring the Gel Pad

The dressing should be changed if the gel pad remains displaced when pressed with a finger. Change Tegaderm™ CHG Dressings every seven days, when the dressing becomes loose or soiled, if the gel pad is saturated, or in cases where there is swelling, visible drainage, or lost visibility.

Removal

1. Using a low and slow removal technique, start removing the dressing from where the catheter or tubing exits the dressing toward the catheter insertion site. Avoid skin trauma by peeling the dressing back, rather than pulling it up from the skin.

2. Apply sterile alcohol swabs or wipes, or sterile solutions (i.e., sterile water or normal saline) between gel pad and skin to facilitate removal of the gel pad dressing.

3. If needed, a medical adhesive solvent can be used to help remove the dressing border. Continue the low and slow removal method until the dressing is completely removed.

When the CHG gel pad is exposed, grasp a corner of the gel pad and the transparent film dressing between thumb and finger.

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